**Monoclonal Mouse  
Anti-Human  
Cytokeratin 19  
Clone RCK108  
Code M0888**

**ENGLISH**

**Intended use**

For in vitro diagnostic use.

Monoclonal Mouse Anti-Human Cytokeratin 19, Clone RCK108, is intended for use in immunohistochemistry (IHC). The antibody labels epithelial cells expressing the cytokeratin 19 protein, and is a useful aid for the classification of epithelial tumors and cholangiocellular carcinoma (1). Differential classification is aided by the results from a panel of antibodies. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist. This antibody is intended to be used after the primary diagnosis of tumor has been made by conventional histopathology using nonimmunologic histochemical stains.

**Summary and explanation**

Cytokeratin 19 (CK 19) belongs to the intermediate filaments, which create a cytoskeleton in almost all cells. In contrast to other intermediate filaments, cytokeratins (CKs) are made up of a highly complex multigene family of polypeptides with molecular masses ranging from 40 to 68 kDa. Twenty distinct CK polypeptides (2, 3) have been revealed in various human epithelial cells (4). They can be divided into an acidic (type I) and a neutral-basic (type II) subfamily. CK 19 belongs to the acidic type of cytokeratins, and is a low molecular mass cytokeratin (40 kDa) typically expressed in simple epithelia, superficial, intermediate and basal cells of transitional epithelium as well as luminal, basal and myoepithelial cells in a few complex epithelia. CK 19 is normally not expressed in stratified squamous epithelia, but may be present in modified squamous epithelium invaded by lymphocytes as well as in basal cells in non-keratinizing stratified squamous epithelium (2, 4, 5).

Refer to *Dako General Instructions for Immunohistochemical Staining* or the detection system instructions of IHC procedures for: Principle of Procedure; Materials Required, Not Supplied; Storage; Specimen Preparation; Staining Procedure; Quality Control; Troubleshooting; Interpretation of Staining; General Limitations.

**Reagent provided**

Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris-HCl, pH 7.2, and containing 15 mmol/L NaNa.

Clone: RCK108 (1). Isotype: IgG1, kappa.

Mouse IgG concentration: See label on vial.

The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot.

**Immunogen**

Cytoskeletal preparation from human bladder cancer cell line T24 (1).

**Specificity**

In Western blotting of cytoskeleton preparations from the T24 and RT4 cell lines as well as the human squamous carcinoma cell line HaCaT, the antibody labels a single band of 40 kDa corresponding to CK 19 (1).

In immunohistochemistry on different human tissues, the antibody labels only epithelia, whereas cells known not to contain CK 19, e.g. hepatocytes, are not labeled (1).

**Precautions**

1. For in vitro diagnostic use.
2. For professional users.
3. This product contains sodium azide (NaN<sub>3</sub>) a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
4. As with any product derived from biological sources, proper handling procedures should be used.
5. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
6. Unused solution should be disposed of according to local, State and Federal regulations.

**Storage**

Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.

**Specimen preparation**

Paraffin sections: The antibody can be used for labeling paraffin-embedded tissue sections fixed in formalin.

Pre-treatment of deparaffinized tissues with heat-induced epitope retrieval (HIER) is required. For heat-induced epitope retrieval, optimal results are obtained with Dako Target Retrieval Solution, pH 9.0, Code S2368, or Dako Target Retrieval Solution, Code S1700. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure.

Frozen sections and cell preparations: The antibody can be used for labeling frozen sections (6). The user must validate the staining procedure.

**Staining procedure**

These are guidelines only. Optimal conditions may vary depending on specimen type and preparation method, and should be validated individually by each laboratory. The performance of this antibody should be established by the user when utilized with other manual staining systems or automated platforms.

Dilution: Monoclonal Mouse Anti-Human Cytokeratin 19, Code M0888, may be used at a dilution range of 1:50-1:100 when applied on formalin-fixed, paraffin-embedded sections of human breast carcinoma and using 20 minutes HIER with Dako Target Retrieval Solution, pH 9.0, Code S2368, and 30 minutes incubation at room temperature with the primary antibody. The recommended negative control is Dako Mouse IgG1, Code X0931, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in Dako Antibody Diluent, Code S0809.

Quality control: Positive and negative control tissues as well as negative control reagent should be run simultaneously using the same protocol as the patient specimens.

Visualization: Dako EnVision+HRP kits, e.g. Code K4005, are recommended for formalin-fixed, paraffin-embedded sections. Follow the procedure enclosed with the selected visualization kit.

**Staining interpretation**

Cells labeled by the antibody display a cytoplasmic staining pattern.

**Performance characteristics**

Normal tissues: The antibody labels basal layer of epidermis (scattered cell labeling), Merkel cells, sebaceous glands, external root sheath of hair follicles, secretory coil and luminal duct cells of eccrine sweat glands, non-keratinizing squamous epithelium of larynx, mesothelium, columnar ciliated epithelial and basal cells of bronchi, alveolar lining cells, cuboidal and basal cells of the parotis gland, ductal, acinar and basal cells of minor salivary glands, basal cells of esophagus, secretory and neuroendocrine cells of the stomach, villus and crypt epithelium of the small intestine, colon epithelium, ductal and acinar epithelium of pancreas, bile duct epithelium, gall bladder epithelium, columnar cells of lobuli and ducti in breast, columnar- and reserve cells of endocervix, basal cells of ectocervix, endometrial glandular epithelium, placental cytotrophoblast,

ductal epithelium, basal cells and seminal vesicles of prostate, rete testis, epithelium of Bowman's capsule, proximal/distal tubuli and Henle's loop, transitional epithelium of the bladder, and follicle epithelium of the thyroid gland (1).

**Abnormal tissues:** In malignant human tissues the antibody labeled: 4/4 basal cell and 3/3 well-differentiated keratinizing squamous cell carcinomas of epidermis, 1/1 keratoacanthoma, 1/3 and 1/2 well-differentiated non-keratinizing squamous cell carcinomas of the lip and tongue respectively, 1/5 moderately and 1/1 poorly differentiated tongue squamous cell carcinomas, 1/1 papilloma and 3/3 inverted papillomas of the nasopharynx, 1/3 well, 3/4 moderately and 1/1 poorly differentiated laryngeal squamous cell carcinomas, 7/7 squamous cell carcinomas, 4/4 adenocarcinomas, 1/1 adenosquamous carcinoma and 5/5 small cell lung carcinomas of the respiratory tract, 5/5 pleomorphic adenomas, 3/3 mucoepidermoid carcinomas and 1/1 acinic cell carcinoma of the parotis gland, 6/6 squamous cell carcinomas of esophagus, 4/4 adenocarcinomas, 2/2 carcinoids and 1/1 lymph node metastasis in the stomach, 1/1 adenocarcinoma of the small intestine, 1/1 carcinoid of the appendix, 11/11 adenocarcinomas, 2/2 lymph node metastasis and 1/1 carcinoid of colon, 7/7 adenocarcinomas of pancreas, 6/6 cholangiocellular carcinomas, 12/12 ductal carcinomas, 4/4 lobular carcinomas, 2/2 lymph node metastasis and 4/4 combined ductal/lobular carcinomas of the breast, 2/2 Paget's disease, 11/11 cystadenocarcinomas, 2/2 clear cell carcinomas and 7/7 endometrioidcarcinomas of the ovary, 1/1 malignant Brenner tumor, 4/4 squamous cell carcinomas, 12/12 adenocarcinomas, 8/8 adenosquamous- and 2/2 clear cell carcinomas of cervix, 13/13 endometrioid carcinomas and 3/3 adenosquamous carcinomas of endometrium, 5/5 well and 1/4 moderately differentiated squamous cell carcinomas of vulva, 8/8 adenocarcinomas of prostate, 3/3 malignant teratomas, 2/2 choriocarcinomas, 6/6 renal cell carcinomas, 9/9 urinary bladder carcinomas, 8/8 thyroid gland carcinomas, 1/1 chromophobic and 1/1 acidophil adenoma of the pituitary gland, 4/4 astrocytomas, 4/4 meningiomas and 5/13 leiomyosarcomas (1). In another study, the antibody labeled 5/8 medullary, 26/28 invasive ductal, 4/4 invasive lobular, 5/5 mucinous and 2/2 tubular malignant breast lesions, as well as 12/14 fibroadenomas, 5/5 Phyllodes tumors, 3/3 tubular adenomas, 1/1 Hamartoma, 3/4 intraduct papillomas, 8/8 gynaecomastias and 5/5 fibrocystic changes of the breast (7). In frozen sections, the antibody labeled 5/5 keratinizing squamous cell carcinomas, 11/11 large-cell non-keratinizing carcinomas, 5/5 endocervical adenocarcinomas, 4/4 endometrioid carcinomas and 3/3 adenosquamous carcinomas of the human cervix (6).

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


1. Kwaspas FHL, Smedts FMM, Broos A, Bulten H, Debie WMH et al. Reproducible and highly sensitive detection of the broad spectrum epithelial marker keratin 19 in routine cancer diagnosis. *Histopathol* 1997;31:503-16.
2. Moll R, Franke WW, Schiller DL, Geiger B, Krepler R. The catalog of human cytokeratins: patterns of expression in normal epithelia, tumors and cultured cells. *Cell* 1982;31:11-24.
3. Moll R, Löwe A, Laufer J, Franke WW. Cytokeratin 20 in human carcinomas. A new histodiagnostic marker detected by monoclonal antibodies. *Am J Pathol* 1992;140:427-47.
4. Moll R. Cytokeratins as markers of differentiation in the diagnosis of epithelial tumors. *Subcell Biochem* 1998; 31: 205-62.
5. Lane EB, Alexander CM. Use of keratin antibodies in tumor diagnosis. In: Osborn M, ed. *Seminars in Cancer Biology. Antibodies in Diagnosis and Therapy*. Philadelphia: Saunders, 1990:165-79.
6. Smedts F, Ramaekers F, Troyanovsky S, Pruszczynski M, Link M, et al. Keratin expression in cervical cancer. *Am J Pathol* 1992;141:497-511.
7. Dalal P, Shousha S. keratin 19 in paraffin sections of medullary carcinoma and other benign and malignant breast lesions. *Mod Pathol* 1995;8:413-6.

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Agilent Technologies Singapore (International) Pte Ltd.