**Monoclonal Mouse  
Anti-Human Epithelial Membrane Antigen**  
Clone E29

**Code M0613**

**ENGLISH**

<b>Intended use</b>	<p>For in vitro diagnostic use.</p> <p>Monoclonal Mouse Anti-Human Epithelial Membrane Antigen, Clone E29, is intended for use in immunohistochemistry (IHC). The antibody labels epithelial cells in a wide variety of tissues and is a useful aid for the classification of neoplastic epithelia (1). Differential classification is aided by the results from a panel of antibodies. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist. This antibody is intended to be used after the primary diagnosis of tumor has been made by conventional histopathology using nonimmunologic histochemical stains.</p>
<b>Summary and explanation</b>	<p>Epithelial membrane antigen (EMA) belongs to a heterogenous population of human milk fat globule (HMFG), proteins. HMFG is a complex secretory product of mammary epithelium and EMA can be recovered from the aqueous phase of skimmed milk following extraction in chloroform and methanol. Besides in milk, these proteins are present in a variety of epithelia of both normal and neoplastic types. EMA is present on the membrane of secretory epithelia. A number of monoclonal and polyclonal antisera have been raised against these molecules and anti-EMA has been extensively studied in a large number of neoplastic conditions in conjunction with other antibodies (2).</p> <p>Refer to <i>Dako General Instructions for Immunohistochemical Staining</i> or the detection system instructions of IHC procedures for: Principle of Procedure; Materials Required, Not Supplied; Storage; Specimen Preparation; Staining Procedure; Quality Control; Troubleshooting; Interpretation of Staining; General Limitations.</p>
<b>Reagent provided</b>	<p>Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris-HCl, pH 7.2, and containing 15 mmol/L NaN<sub>3</sub>.</p> <p><u>Clone:</u> E29 (1). <u>Isotype:</u> IgG2a, kappa.</p> <p><u>Mouse IgG concentration:</u> see label on vial.</p> <p>The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot.</p>
<b>Immunogen</b>	Human milk fat globule membrane preparation (1, 3).
<b>Specificity</b>	In Western blotting of the immunogen the antibody labels bands of 265-400 kDa (1).
<b>Precautions</b>	<ol style="list-style-type: none"> <li>1. For in vitro diagnostic use.</li> <li>2. For professional users.</li> <li>3. This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.</li> <li>4. As with any product derived from biological sources, proper handling procedures should be used.</li> <li>5. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.</li> <li>6. Unused solution should be disposed of according to local, State and Federal regulations.</li> </ol>
<b>Storage</b>	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Support.
<b>Specimen preparation</b>	<p><u>Paraffin sections:</u> The antibody can be used for labeling paraffin-embedded tissue sections fixed in formalin, Bouin's fixative, B5 fixative or Zenker's solution (5). Pre-treatment of deparaffinized tissues with proteinase K or heat-induced epitope retrieval is recommended. For heat-induced epitope retrieval of tissues fixed in formalin, optimal results are obtained with Dako Target Retrieval Solution, Code S1700, 10 mmol/L citrate buffer, pH 6.0, or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure.</p> <p><u>Frozen sections and cell preparations:</u> The antibody can be used for labeling acetone-fixed frozen sections (1). The user must validate the staining procedure.</p>
<b>Staining procedure</b>	<p>These are guidelines only. Optimal conditions may vary depending on specimen type and preparation method, and should be validated individually by each laboratory. The performance of this antibody should be established by the user when utilized with other manual staining systems or automated platforms.</p> <p><u>Dilution:</u> Monoclonal Mouse Anti-Human Epithelial Membrane Antigen, Code M0613, may be used at a dilution range of 1:50-1:100 when applied on formalin-fixed, paraffin-embedded sections of human tonsil and using 20 minutes heat-induced epitope retrieval in Dako Target Retrieval solution, Code S1700, and 30 minutes incubation at room temperature with the primary antibody. The recommended negative control is Dako Mouse IgG2a, Code X0943, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in Dako Antibody Diluent, Code S0809.</p> <p><u>Quality control:</u> Positive and negative control tissues as well as negative control reagent should be run simultaneously using the same protocol as the patient specimens..</p> <p><u>Visualization:</u> Dako EnVision+/HRP kits, e.g. Code K4005, are recommended. Follow the procedure enclosed with the selected visualization kit.</p>
<b>Product-specific limitations</b>	The antibody labels plasma cells, and EMA appears to be common in plasma cell neoplasms, but is also encountered occasionally among other types of lymphoma. However, it should be emphasized that such cases usually are clearly identifiable on purely morphological grounds as being of lymphoid origin (1).
<b>Staining interpretation</b>	In normal breast and other secretory epithelia, labeling is predominantly localized to apical luminal membranes. In neoplasms, cytoplasmic and apical luminal membrane staining are the most common patterns of immunoreactivity with peripheral membrane staining or other patterns also occurring (4).

**Performance characteristics** Normal tissues: The antibody labels epithelial cells in a wide variety of tissues and mesothelial cells. Included are sweat ducts and sebaceous glands of the skin, epithelium of the gastrointestinal tract, acini and ducts of the breast, exocrine cells in pancreas, bladder epithelium, distal tubules of the kidney, cervix, endometrium, respiratory epithelium, thyroid and bile ducts. No labeling was observed in epidermis, endocrine cells in pancreas, glomeruli and proximal tubules of the kidney, central nervous system, peripheral nervous system, connective tissue, hepatocytes, and lymphoid tissue, except for occasional plasma cells (1).

Abnormal tissues: The antibody labels a wide variety of neoplastic epithelia, and also neoplastic mesothelial cells (1). In a large study of 2081 epithelial, mesenchymal and hematopoietic neoplasms, the antibody provided a 98.6% specificity (2). It was shown that the antibody labeled 105/354 cases of soft tissue and intracranial tumors, with most labeled cases among synovial sarcomas, spindle cell malignant mesotheliomas, epithelioid sarcomas, chordomas and choroid plexus tumors, 38/169 cases of small round cell tumors and small cell sarcomas, 13/158 cases of germ cell neoplasms, 23/23 cases of spindle cell "sarcomatoid" carcinomas, 815/918 cases of other epithelial malignancies, including squamous, Merkel cell, breast, gastric adeno-, colonic adeno-, pancreatic, salivary gland, bladder, uterine, ovarian, vaginal, pulmonary, prostatic, thyroid, thymic, hepatic, renal cell, and nasopharyngeal carcinomas (2). In lymphomas, the antibody labeled 10/22 cases of null or T-cell type CD30+ ALCL of small cell, monomorphic, or pleomorphic subtypes (5), 34/35 cases of p80/ALK+, 1/6 cases of CD56/57+ T/NK-cell, 2/7 cases of EBV+ cytotoxic large T-cell, 2/8 cases of low-grade cytotoxic T-cell and 6/10 cases of cytotoxic Hodgkin's-like lymphomas (6). No labeling was observed in 18/18 cases of basal cell carcinoma, 12/12 cases of hepatocellular carcinoma, 43/43 cases of malignant melanoma and 69/69 cases of endocrine neoplasms, with the exception of poorly differentiated lesions in 6 cases of islet cell tumors (2).

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


1. Cordell J, Richardson TC, Pulford KAF, Ghosh AK, Gatter KC, Heyderman E, et al. Production of monoclonal antibodies against human epithelial membrane antigen for use in diagnostic immunocytochemistry. Br J Cancer 1985;52:347-54.
2. Swanson PE. Monoclonal antibodies to human milk fat globule proteins. In: Wick MR, Siegal GP, editors. Monoclonal antibodies in diagnostic immunohistochemistry. New York – Basel: Marcell Dekker Inc; 1988. p. 227-83.
3. Heyderman E, Strudley I, Powell G, Richardson TC, Cordell JL, Mason DY. A new monoclonal antibody to epithelial membrane antigen (EMA) – E29. A comparison of its immunocytochemical reactivity with polyclonal anti-EMA antibodies and with another monoclonal antibody, HMF2-2. Br J Cancer 1985;52:355-61.
4. Pinkus GS, Kurtin PJ. Epithelial membrane antigen – a diagnostic discriminant in surgical pathology. Immunohistochemical profile in epithelial, mesenchymal, and hematopoietic neoplasms using paraffin sections and monoclonal antibodies. Hum Pathol 1985;16:929-40.
5. Felgar RE, Salhany KE, Macon WR, Pietra GG, Kinney MC. The expression of TIA-1+ cytolytic-type granules and other cytolytic lymphocyte-associated markers in CD30+ anaplastic large cell lymphomas (ALCL): correlation with morphology, immunophenotype, ultrastructure, and clinical features. Hum Pathol 1999;30:228-36.
6. Kagami Y, Suzuki R, Taji H, Yatabe Y, Takeuchi T, Maeda S, et al. Nodal cytotoxic lymphoma spectrum. A clinicopathologic study of 66 patients. Am J Surg Pathol 1999;23:1184-1200.

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Agilent Technologies Singapore (International) Pte Ltd.