**Polyclonal Rabbit
Anti-Human
Chorionic Gonadotropin**

Code A0231

ENGLISH

Intended use	<p>For in vitro diagnostic use.</p> <p>Polyclonal Rabbit Anti-Human Chorionic Gonadotropin (hCG), Code A0231, is intended for use in immunohistochemistry (IHC). The antibody labels hCG-containing cells (1) and may be used for the demonstration of trophoblastic elements. Results aid in the classification of germ cell tumors (2). Differential classification is aided by the results from a panel of antibodies. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist. This antibody is intended to be used after the primary diagnosis of tumor has been made by conventional histopathology using nonimmunologic histochemical stains.</p>
Summary and explanation	<p>hCG is a 40 kDa glycoprotein composed of two subunits, α and β. The α subunit is common for hCG, LH, follicle stimulating hormone (FSH), and thyroid stimulating hormone (TSH), while the β subunit is unique for hCG although its amino acid sequence partially is identical with that of the β subunit of LH (3, 4). hCG is the most important marker of gestational trophoblastic cells, being present in syncytiotrophoblastic cells and cells of the intermediate trophoblast but absent in cytotrophoblasts. hCG expression has also been demonstrated in some non-trophoblastic tumors, e.g. in some patients with hepatocellular carcinoma (4).</p> <p>Refer to Dako General Instructions for Immunohistochemical Staining or the detection system instructions of IHC procedures for: Principle of Procedure; Materials Required, Not Supplied; Storage; Specimen Preparation; Staining Procedure; Quality Control; Troubleshooting; Interpretation of Staining; General Limitations.</p>
Reagent provided	<p>Purified immunoglobulin fraction of rabbit antiserum provided in liquid form. In 0.1 mol/L NaCl, 15 mmol/L Na₂S₂O₃.</p> <p><u>Protein concentration:</u> See label on vial.</p> <p>The titre variation between different lots of A0231 as determined by single radial immunodiffusion is less than 10%. This is achieved by adjusting the titre of each individual lot to match the titre of an antibody reference preparation kept at -80 °C.</p>
Immunogen	<p>The isolated beta-chain of hCG.</p>
Specificity	<p>The antibody reacts with hCG, and like most other polyclonal antibodies to beta-hCG it also shows a slight cross-reaction (about 4%) with LH (3). Traces of contaminating antibodies to human plasma proteins have been removed by solid-phase absorption.</p> <p>The specificity of the antibody has been tested as follows:</p> <p><u>Crossed immunoelectrophoresis:</u> Only the hCG precipitation arch appears when using 12.5 μL A0231 per cm² gel area against 10 μL of concentrated urine from pregnant women (15th-16th week of gestation). No precipitation is seen when 12.5 μL antibody per cm² gel area is tested against 10 μL of normal human plasma. Staining: Coomassie Brilliant Blue.</p>
Precautions	<ol style="list-style-type: none"> 1. For in vitro diagnostic use. 2. For professional users. 3. This product contains sodium azide (Na₃N), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 4. As with any product derived from biological sources, proper handling procedures should be used. 5. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. 6. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	<p>Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, relevant controls should be run simultaneously with patient specimens. If unexpected results are observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.</p>
Specimen preparation	<p>The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin. Pre-treatment of deparaffinized tissues with heat-induced epitope retrieval is required. Optimal results are obtained with Dako Target Retrieval Solution, Code S1700, or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. Less optimal results are obtained with 10 mmol/L citrate buffer, pH 6.0 and pre-treatment of tissues with proteinase K. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure.</p>
Staining procedure	<p>These are guidelines only. Optimal conditions may vary depending on specimen type and preparation method, and should be validated individually by each laboratory. The performance of this antibody should be established by the user when utilized with other manual staining systems or automated platforms.</p> <p><u>Dilution:</u> Polyclonal Rabbit Anti-Human Chorionic Gonadotropin (hCG), Code A0231, may be used at a dilution range of 1:300-1:600 when applied on formalin-fixed, paraffin-embedded sections of human placenta and using 20 minutes heat-induced epitope retrieval in Dako Target Retrieval Solution, Code S1700, and 30 minutes incubation at room temperature with the primary antibody. The recommended negative control is Dako Rabbit Immunoglobulin Fraction (Solid-Phase Absorbed), Code X0936, diluted to the same protein concentration as the primary antibody. Unless the stability in the actual test system has been established, it is recommended to dilute the product immediately before use or dilute in Dako Antibody Diluent, Code S0809.</p> <p><u>Quality control:</u> Positive and negative control tissues as well as negative control reagent should be run simultaneously using the same protocol as the patient specimens.</p> <p><u>Visualization:</u> Dako EnVision+/HRP kits, e.g. Code K4009, are recommended. Follow the procedure enclosed with the selected visualization kit.</p>
Staining interpretation	<p>Cells labelled by the antibody display a cytoplasmic staining pattern (1).</p>

Performance characteristics

Normal tissues: In first trimester placenta, the antibody labels plentiful syncytiotrophoblast cytoplasmic granules, some undergoing exocytosis (1).

Abnormal tissues: 189 orchidectomy specimens with germ cell tumors (95 pure seminomas and 94 non-seminomas) were studied with the antibody. hCG was demonstrated in all choriocarcinoma components in the syncytiotrophoblasts and syncytiotrophoblast-like cells in 8% of seminomas and 30% of non-seminomas. Carcinoma-in-situ present in seminiferous tubules adjacent to the tumors was labeled for hCG in a few cases (2).

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
1. Morrish DW, Marusyk H, Siy O. Demonstration of specific secretory granules for human chorionic gonadotropin in placenta. J Histochem Cytochem 1987;35:93-101.
2. Jacobsen GK, Jacobsen M. Alpha-fetoprotein (AFP) and human chorionic gonadotropin (HCG) in testicular germ cell tumours. A prospective immunohistochemical study. Acta Path Microbiol Immunol Scand Sect A 1983;91:165-76.
3. Mehta HC, MacDonald DJ. A sensitive enzyme immunoassay specific for human chorionic gonadotrophin. Clin Chim Acta 1982;121:245-50.
4. Leong AS-Y, Cooper K, Leong FJW-M. Manual of diagnostic antibodies for immunohistology; London: Oxford University Press; 1999. p. 25-6.

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Agilent Technologies Singapore (International) Pte Ltd.