**Monoclonal Mouse  
Anti-Vimentin**  
Clone V9

**Code M0725**

**ENGLISH**

<b>Intended use</b>	<p>For in vitro diagnostic use.</p> <p>Monoclonal Mouse Anti-Vimentin, Clone V9, is intended for use in immunohistochemistry (IHC). The antibody labels primarily cells of mesenchymal origin in normal and neoplastic tissues, and is a useful aid for classification of tumors of mesenchymal origin (1). Differential classification is aided by the results from a panel of antibodies. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist. This antibody is intended to be used after the primary diagnosis of tumor has been made by conventional histopathology using nonimmunologic histochemical stains.</p>
<b>Summary and explanation</b>	<p>Vimentin is a 57 kDa intermediate filament (IF) protein, which form part of the cytoskeleton of vertebrate cells. Among the IFs, vimentin belongs to class III, showing a high degree of specificity for cells of mesenchymal origin. The cell type specificity, displayed by each of the IF subtypes, was initially thought to be retained in malignant cells as well as their normal counterpart, which made IFs important as markers in histogenesis. The coexpression of intermediate filaments, particularly vimentin and cytokeratin, has been demonstrated in a variety of normal cells/tissues and in neoplastic lesions, necessitating the use of a panel of antibodies in tumor classification (2).</p> <p>Refer to <i>Dako General Instructions for Immunohistochemical Staining</i> or the detection system instructions of IHC procedures for: Principle of Procedure; Materials Required, Not Supplied; Storage; Specimen Preparation; Staining Procedure; Quality Control; Troubleshooting; Interpretation of Staining; General Limitations.</p>
<b>Reagent provided</b>	<p>Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris-HCl, pH 7.2, and containing 15 mmol/L Na<sub>2</sub>S<sub>2</sub>O<sub>3</sub>.</p> <p><u>Clone:</u> V9 (3). <u>Isotype:</u> IgG1, kappa.</p> <p><u>Mouse IgG concentration:</u> see label on vial.</p> <p>The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot.</p>
<b>Immunogen</b>	<p>Purified vimentin from porcine eye lens (3).</p>
<b>Specificity</b>	<p>In Western blotting of purified porcine vimentin, the antibody labels a single band of 57 kDa corresponding to vimentin. When applying whole cell extracts of cell lines expressing vimentin plus glial fibrillary acidic protein (GFAP), and vimentin plus desmin, respectively, the antibody labels specifically the 57 kDa vimentin band. As directly shown by these experiments, the antibody does not react with the two IF proteins most closely related to vimentin, i.e. desmin and GFAP (3).</p> <p>In IHC, the antibody labels the vimentin-positive human cell lines IMR90, RD, glioma and HeLa (3).</p>
<b>Precautions</b>	<ol style="list-style-type: none"> <li>1. For in vitro diagnostic use.</li> <li>2. For professional users.</li> <li>3. This product contains sodium azide (Na<sub>2</sub>N<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.</li> <li>4. As with any product derived from biological sources, proper handling procedures should be used.</li> <li>5. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.</li> <li>6. Unused solution should be disposed of according to local, State and Federal regulations.</li> </ol>
<b>Storage</b>	<p>Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.</p>
<b>Specimen preparation</b>	<p><u>Paraffin sections:</u> The antibody can be used for labeling paraffin-embedded tissue sections fixed in formalin. Pre-treatment of deparaffinized tissues with heat-induced epitope retrieval is required. Optimal results are obtained with 10 mmol/L citrate buffer, pH 6.0. Less optimal results are obtained with 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. However, Dako Target Retrieval Solution, Code S1700, and pre-treatment of tissues with proteinase K were found destructive of the epitope. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure.</p> <p><u>Frozen sections and cell preparations:</u> The antibody can be used for labeling frozen sections or fixed cell smears (3). The user must validate the staining procedure.</p>
<b>Staining procedure</b>	<p>These are guidelines only. Optimal conditions may vary depending on specimen type and preparation method, and should be validated individually by each laboratory. The performance of this antibody should be established by the user when utilized with other manual staining systems or automated platforms.</p> <p><u>Dilution:</u> Monoclonal Mouse Anti-Vimentin, Code M0725, may be used at a dilution range of 1:50-1:100 when applied on formalin-fixed, paraffin-embedded sections of human tonsil and using 20 minutes heat-induced epitope retrieval in 10 mmol/L citrate buffer, pH 6.0 and 30 minutes incubation at room temperature with the primary antibody. The recommended negative control is Dako Mouse IgG1, Code X0931, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in Dako Antibody Diluent, Code S0809.</p> <p><u>Visualization:</u> Dako EnVision+HRP kits, e.g. Code K4005, are recommended. Follow the procedure enclosed with the selected visualization kit.</p> <p><u>Quality control:</u> Positive and negative control tissues as well as negative control reagent should be run simultaneously using the same protocol as the patient specimens.</p>
<b>Staining interpretation</b>	<p>Cells labeled by the antibody display a cytoplasmic staining pattern.</p>

<b>Performance characteristics</b>	<p><b>Normal tissues:</b> In general, most human mesenchymal cells are labeled by the antibody, including fibrocytes, lipocytes, smooth muscle cells, vascular endothelial cells, astrocytes, peripheral nerve (Schwann) cells, macrophages (including Kupffer cells), as well as myoepithelial cells of sweat and salivary glands and of breast, which are all labeled strongly. Also labeled, with variable intensity and distribution, are the follicular cells of the thyroid, adrenal cortex, renal distal tubules, and mesangial and endothelial cells of the renal glomerulus, as well as pancreatic acinar cells (1, 3). In the human eye, the antibody labels the pigmented posterior and the anterior epithelia of the human iris, including the muscle portion (dilator pupillae) of the anterior epithelium, as well as the nonpigmented and pigmented ciliary epithelia (4). In the ciliary epithelia, vimentin was coexpressed with cytokeratin (4). Skeletal and cardiac muscle cells, epidermal, squamous, urothelial, colonic and gastric mucosal, and glial cells, as well as neurons are consistently not labeled with the antibody (1, 3).</p> <p><b>Abnormal tissues:</b> The antibody labeled 17/20 sarcomas, 16/18 melanomas, 4/4 meningiomas, and 3/3 schwannomas, and was the sole intermediate filament present in these tumors. In addition, variable percentages (10-57%) of carcinomas, neuroendocrine carcinomas, neuroblastomas, thymomas and mesotheliomas were labeled with the antibody. With the exception of the neuroblastomas, cytokeratin was coexpressed with vimentin in these tumors. Among adenocarcinomas, more than 50% of papillary carcinomas of the thyroid as well as renal, endometrial, ovarian and lung carcinomas were labeled by the antibody and coexpressed keratins and vimentin (1).</p>
------------------------------------	--

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Agilent Technologies Singapore (International) Pte Ltd.